

SM-500S

12.1inchTFT



Picture:

- Big font display interface
- Trend graph display interface
- OxyCRG display interface
- ECG, NIBP, SPO2, TEMP interface
- Standard display interface
- Built-in rechargeable battery and recorder
- ECG full lead display interface

Feature:

- Display with 9 waveform
- 720 hours trend data
- 1000 NIBP data storage
- 200 alarm event data storage
- Sound, light, human voice alarm
- 12 hours waveform review
- Wire and wireless network
- Nurse call
- analog output
- Support software upgrade online
- Arrhythmia analysis, ST segment analysis
- Support touch screen

Standard:

ECG, RESP, NIBP, SpO2, TEMP

Optional:

Recorder, IBP, EtCO2, touch screen



Sonostar Technologies Co., Limited

Address: #31, Shuiwan road, Jida, Zhuhai, Guangdong, China
Tel: +86 13168355297 MSN: sonostar@163.com Skype: sonostar-ut
Web: www.sonostar.net Mail: sonostar@sonostar.net





SM-500S



- 12.1 inch color high resolution TFT display
- Portability, suitable for at patient's bedside, ICU monitor
- Unique blood pressure measurement technique
- Anti high-frequency electro-tome and defibrillation, digital filter
- Built-in rechargeable battery, DC/AC

